



EZ BreathSure

Visible Display

- ▶ Large-capacity battery, long working time for 40 hours;
- ▶ Built-in Lithium battery, Type-C interface charging;
- ▶ LED digital tube displays respiratory rate and alarm;
- ▶ Available to measure respiratory rate;
- ▶ Apnea alarm, Exceed-limit alarm;
- ▶ Apnea with Audible alarm;
- ▶ Respiratory rhythm sound;
- ▶ Automatic wake up breathing, automatic shutdown without breathing detection.



Product Specifications

Physical Specifications	
Size	• 68.7 mm (h) x 62.85mm (w) x 24.2mm (d) • (2.7 in (h) x 2.47 in (w) x 0.95in (d))
Weight	30g
Power Supply	
External Power	5V DC, 0.4A, Type-C
Battery Specifications	
Battery capacity	3.7V/350mAh
Operating Time	Operating Time \geq 40hours

Product Configuration	
Respiratory Rate, Apnea Alarm, Li-On rechargeable Patient type: Cat, Dog, Horse	
Environment Specification	
Temperature	0°C ~ 40°C
Humidity	15% ~ 80%, non-condensing

Environment Specification	
Accuracy	\pm 2rpm
Measurement Range	0~120rpm
Alarm Limits	>50rpm or <5rpm
Apnea Alarm	20s
Volume	>70dB

Display	
Display	3-segment white LED

Tootoo Meditech Co., Ltd.

3A, Phase II, Smart Home, Baolong Street, Longgang District, Shenzhen, Guangdong, China

+86-755-28306385

www.newtootoo.com

info@newtootoo.com