



## DTS-180 Animal Dental System

Simplify Animal Oral Treatments Efficiently

Provides a variety of functions for veterinarians' clinic needs, grinding, polishing, cutting, cleaning etc



### Stable Compressor

- Wrapped aluminum alloy heat sink.
- Independent gas circuit design and easy operation
- Explosion-proof tube over pressure nozzles
- Air source : Air pressure 0.8Mpa



### Standard Hand Pieces

#### AJ-T1

- Single spray
- Ceramic bearing
- Push button type
- Open cartridge
- 4 holes turbine
- 134 C autoclavable
- 320000-4000rpm high-frequency scaling handle

#### AJ-S1

- A bend low-speed hand piece, straight head in option
- 134 C autoclavable
- Rotate speed: 20000rpm
- Air Pressure: 300-350Kpa
- Atomize Air/water pressure: 200K
- CA Burs : 2.35mm.