




Veterinary Weight Scale TTPS-3624



PRODUCT DETAILS



EASY WEIGHING

Special made for animals, large platform and easy to weigh



LED SCREEN

Accurate measurement, easy reading



ANTI-SLIP MAT

Considerate design, prevent falling from smooth surface



ADJUSTABLE FEET

Four adjustable feet provide excellent stability

WALK-ON SCALE

- Simple shape and excellent detail appearance
- Made of 304 stainless steel, easy to clean
- LED display with cable (Bluetooth is optional)
- Anti-slip mat

PRODUCT PARAMETER

Product	walk-on scale
Size	60*100cm
Scaling Range	0.5kg-300kg
Accuracy	0.05-0.1KG
Power	Operate with both AC adapter and battery

Tootoo Meditech Co., Ltd.

Address: 3A, Phase II, Smart Home, Baolong Street, Longgang District, Shenzhen, Guangdong, China

Tel: +86-755-28306385

Email: info@newtootoo.com

Web: www.newtootoo.com