

Digital Video Monitor



***Simultaneous Dual-Scope
Display for Enhanced
Diagnostic Capabilities***

Hunan Vathin Medical Instrument Co.,Ltd.

1/F, Building 12, Innovation and Entrepreneurship Service Center No. 9
Chuanqi West Road, Jiuhua Economic Development Zone,
411100 Xiangtan, Hunan, China

 www.vathin.com  marketing@vathin.com

Application:

Pulmonology, Urology, ENT...

Product features:

Display: Large screen size for easier image examination and recognition.

Strong Compatibility: Compatible with Vathin endoscopes of various specifications.

Connectivity: Seamless integration with hospital PACS system.

Ample Power Supply: Equipped with a high-capacity rechargeable lithium battery.




Image Quality: Built-in image adjustment features with multiple department-specific image modes.

Remote Control: for video transmission via HDMI/DVI-D to PACS.

Enhanced Clarity: Super-resolution imaging.

Versatile Use: Compatible with various other energy devices.

Specifications:

Model	DVM-D1
Orientation:	Landscape, Portrait
Display resolution:	1920*1080 (full HD 60Hz)
Playback:	Picture(connect to WiFi for direct printing) & Video
Dimension WxHxT:	382-mm x 238mm x 44.5mm
Weight	3 kg
Storage Capacity	64G
Display Type	15.6" FHD Color TFT-LCD
Image Adjustment	Color, Contrast, Sharpness, Brightness, and Gamma
Battery Type	10.8V 9000mAh Li-ion
Battery Operation	More than 4 hours
File transfer	USB, DICOM
USB Connection	2 x USB 3.0 Type A
Video output	HDMI & SDI(3G), 1920X1080P60