

Digital Video Monitor – B1



DVM

Specification

Video Record Duration more than 6 hours

Playback Picture & Video

Electrical Power Requirement

100V-240V AC; 47-63HZ; 1.62-0.72A

Dimension WxHxT 310mmx 231mm x 39mm

Weight 2.14kg **Memory** 64G

Function



- ⊙ Max Resolution: (16:10) -1280 x 800 Pixels
- ⊙ Display Type: 12.1" Colour TFT LCD
- ⊙ Image Adjustment: Sharpness and Brightness, Gamma parameters, RGB Saturation
- ⊙ Battery Type: 10.8V 9000mAh Li-ion
- ⊙ Battery Operation: More than 6 hours
- ⊙ File transfer: USB Connection: 2 x USB 3.0 Type A
- ⊙ Video out: HDMI (1920 x 1080p, 60 fps)
- ⊙ PACS Connectivity: Yes (Wi-Fi&LAN)