

BMO210/310 Vet

Vital Signs Monitor Intuitive. Adaptable. Reliable.

The BMO210/310 Vet can be used on small, medium, or large animals. An applicable monitor for veterinary surgical suites and research facilities. Its 12-inch screen has a straightforward and predictable interface, enabling you to easily perform functions and allows for easy-to-view vital signs monitoring at a glance. For your convenience, the intuitive monitors can be controlled via touchscreen, knob or hard keys, and patient vital signs data download via Wi-Fi or Ethernet for each case. Integrated printers are available on monitors including basic and advanced parameters. As your needs change, you can easily upgrade your BMO210/310 Vet monitor in your practice.



Features

Large BMO210 12.1"/BMO310 15" display provides a brighter screen with better contrast and better visual effects;

Large font display to offer clear view of all vital signs;

Large storage function for NIBP measurement, alarm event, trend diagram/table review, etc;

Enhanced storage compared to other leading monitor brands: 48 hours waveform data 1000 hours trend diagram review 1000 hours trend table review 120 minutes dynamic trend 1000 NIBP measurements 200 alarm events.

Multiple Configurations Available:

Standard Parameters: ECG, Pulse Oximetry, Temperature, NIBP, Respiration, Optional Parameters: Invasive Blood Pressures, EtCO₂ (Mainstream or Sidestream), WIFI, Suntech NIBP/IBP, Nellcor Spo₂, Masimo Spo₂, Touchscreen, Thermal Recorder; Comes with standard veterinary accessories;
Versatile mounting options : bedside hook, wall-mount moving cart



BMO210A

12.1" inch Monitor



BMO210B

Standard with Sidestream ECO₂



BMO310A

15 inch Monitor



BMV Medical Group Co., Ltd

Manufacturer of specialty medical devices across
medical and veterinary fields since 2008.

Shenzhen Qianhai headquarters
Add: 4/F, Yinjin Building, Block 71, Baoan Centre District,
Shenzhen, China
Tel: +86 75526564580
W: www.bmv.cc www.bmv-vet.com
Hong Kong Branch
Address: 5 / F, Lianwei Commercial Building, 283 Queen's
Road Central, Hong Kong
Tel: (00852)28510669
Email: export03@bmvc.com

