




TooToo MediTech Co., Ltd.

TT-FC100

Vet Fundus Cameras

Handheld camera offers easy-to-use for veterinary retinal imaging and documentation.



Address: 3A, Phase II, Smart Home, Baolong Street,
Longgang District, Shenzhen, Guangdong, China
Tel: +86-755-28306385 Email: info@newtootoo.com
Web: www.newtootoo.com



42° wide-angle viewing



Built-in battery supports 4-6 hours working



Non-mydriatic



WiFi transmission for easy data sharing



HD imaging with CMOS sensor



Portable and handheld



Specifications

Size	19cm*10cm*21cm
Net Weight	400g (Typical)
Field of View	42 degree
Minumum Pupil	2.5mm
Image Resolution	1920*1080(1080P)
Refractive Compensation	-20D~+20D
Focus Mode	Manual
Flash Mode	White LED/ Near Infrared LED
LCD Screen	3.5 inch TFT Colored Screen
Picture Format	JPEG(Photograph)
Interface	Micro USB
Image Transmission	Image Transmission: Micro USB , Wifi Transmission
Storage	Micro SD card. Maximum to 32G
Power	Li-on rechargeable Battery 3.6V/3400mAh