

Features

- Strong plastic case with metal lid
- Accurately measuring the specific value of blood cell
- Separating micro-amount of blood or solution
- With reader
- LED display speed and timer
- Error code display, speed up and down adjustment
- With mechanical door lock and safety switch
- With brake and automatic alarm

Technical Specifications

Model	HC-12C
Max RPM	12000RPM
Max RCF	15300xg
Timer	0-99min
Rotor capacity	Φ1.5x75mmx24 capillary tubes
Voltage	110V 60Hz /220V 50Hz
Dimensions	280*310*265mm
Package	340*350*330mm
G.W./N.W.	10/9KG

