



深圳市谊业仪器有限公司

T00T00 MEDITECH CO., LTD.

Address: RM 607, HAIHANG CITY BUSINESS CENTER, LONGGANG, SHENZHEN, CHINA.

Web.: www.newtootoo.com

Tel: +86-755-28306385

Email: info@newtootoo.com

MK-1206A 6 Channel ECG Specifaicon

Lead	12 leads,Standard or Cabrera, 12 leads simultaneously acquisition
Input circuit	Floating,Defibrillation protection
LCD display	800*480 7 inch touchscreen LCD display ECG waveform, operation conditions,Time,Heart-rate,etc
Safety	IEC class I, type CF
Sampling rate	1000 samples/Sec per channel
Input dynamics	±500mV
Operation Mode	automatic,manual,single channel 60s compressed record,Rhythm analysis (R-R histogram and trend chart)
Filter	AC(50Hz or 60Hz, -20dB)
	EMG (25Hz/35Hz/45Hz/75Hz/100Hz,-3dB)
	Drift (0.5Hz,-3dB)
CMRR	>100dB(with AC filter)
Input circuit current	≤0.1μA
Input Impedance	≥50Mohm
Time constant	≥ 3.2 Sec
Frequency Response	0.05--150Hz,-3dB
Noise Level	≤15μVp-p
Sensitivity	1.25; 2.5; 5; 10;20; 40mm/mV±2%
Calibration Voltage	1mV±2%
Patient Current Leakage	<10uA
Protection	Protect from defibrillation and pace-making
Printing system	Thermal printer; ≥8dot/mm(Vertical), ≥40 dot/mm (Horizontal),printing speed at 25mm/s
Printing Speed	5; 6.25; 10; 12.5; 25; 50mm/s±3%
Printing Output	3channel; 3+1 channel; 6channel; 6+1 channel
Paper size	110mm*20m roll type or Z-fold thermal paper
Patient Cable	standard 10 core
Report storage	200 reports in built-in memory, support external SD card
	report presentation,review and re-printing.
Indicator	Lead-off detectoin
Power Supply	12.5V 1800mA Li-poly rechargeable battery pack (DC)
	220V±10%, 50Hz±2%, 75VA (AC)
Fuse Specification	2Φ5x20mm 3.15A 250V
Machine Dimension	285(W)*200(D)*55mm(H)
Machine weight	2.25 kg (machine only)
Gross weight	3.9 kg (with accessories)
Carton size	41*32*24cm(1pc), 67*44*51cm(4pcs/carton)

MK-1206A

7 inch touchscreen ECG machine
(6 channel)

For Veterinary Use Only

